

Lactoflow

Generic Name	Fennel Extract, Fenugreek Extract, Cumin powder, Cumin Extract, Dill Powder, Aniseed Extract
Brand Slogan	Induces Lactation Naturally
Composition	Each Lactoflow capsule contains Fennel Extract 83.33 mg, Fenugreek Extract 83.33 mg, Cumin powder 83.33 mg, Cumin Extract 83.33 mg, Dill Powder 16.67 mg and Aniseed Extract 16.67 mg
Description	Breast milk is considered to be an ideal and the optimal food source for newborns through one year of age. It contains everything the baby needs for the first six months of life. The composition of breast milk even changes according to the baby's changing needs. Breastfeeding promotes healthy weight gain and helps prevent childhood obesity. Lactation failure due to inadequate breast feeding or insufficient breast milk production can cause several problems in the initial days of life. The mortality and morbidity risks of infants who receive other nutrition are four times more than infants who are breast fed. According to the World Health Organization (WHO) and UNICEF 1.5 million babies die every year because they were not breastfed. Many more millions of babies suffer from infectious diseases and malnutrition, never reaching their full potential because they were bottle-fed. Breast feeding is the best nutrition for infants. Most of the postnatal mothers <u>face</u> the problem of insufficient breast milk which affects the health and well-being of an infant.
Indication	Lactoflow capsule is indicated for increased breast milk production in lactating mothers.
Dosage & Administration	1 capsule 3 times daily after meal or as advised by physician.
Side Effects	No significant side-effects have been found so far in the prescribed dose. Some side-effects like nausea, vomiting, diarrhea, gastric irritation etc may be happened.
Contraindication	The herbs of Lactoflow capsule is contraindicated in patients who are hypersensitive to any of its ingredients.

Use in during pregnancy and lactation

Lactoflow capsule should not be used during pregnancy. Because Ingredients used in the **Lactoflow** capsule can induce unintentional labor.

Storage

Store in a cool and dry place below 30° C temperature, Keep away from light, heat and moisture.

Presentation

Each box contains 30's Capsules in Alu-Alu blister pack.